

Biotin-PEG4-alkyne

Chemical Properties

CAS No. : 1458576-00-5


Formula: C₂₁H₃₅N₃O₆S

Molecular Weight: 457.58

Store at low temperature

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Biotin-PEG4-alkyne is a biotin-labeled, polyethylene glycol-based (PEG-based) PROTAC linker utilized in PROTAC synthesis.
Targets(IC50)	PROTAC Linker
In vitro	PROTACs contain two different ligands connected by a linker; one of them is a ligand for an E3 ubiquitin ligase and the other one is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Solubility Information

Solubility	DMSO: 250 mg/mL (546.35 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 4 mg/mL (8.74 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1854 mL	10.9271 mL	21.8541 mL
5 mM	0.4371 mL	2.1854 mL	4.3708 mL
10 mM	0.2185 mL	1.0927 mL	2.1854 mL
50 mM	0.0437 mL	0.2185 mL	0.4371 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Gadd MS, et al. Structural basis of PROTAC cooperative recognition for selective protein degradation. Nat Chem Biol. 2017 May;13(5):514-521.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481