


Notoginsenoside FP1

Chemical Properties

CAS No. :	1004988-73-1
Formula:	C47H80O18
Molecular Weight:	933.14
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Notoginsenoside FP1 is a natural product that can be used in related research in the field of life sciences. Its product number is TN8463.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0717 mL	5.3583 mL	10.7165 mL
5 mM	0.2143 mL	1.0717 mL	2.1433 mL
10 mM	0.1072 mL	0.5358 mL	1.0717 mL
50 mM	0.0214 mL	0.1072 mL	0.2143 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481