

Isochiisanoside methyl ester

Chemical Properties


CAS No. : 105273-99-2

Formula: C49H78O20

Molecular Weight: 987.14

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Isochiisanoside methyl ester is a natural product that can be used in related research in the field of life sciences. Its product number is TN8473.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.013 mL	5.0651 mL	10.1303 mL
5 mM	0.2026 mL	1.013 mL	2.0261 mL
10 mM	0.1013 mL	0.5065 mL	1.013 mL
50 mM	0.0203 mL	0.1013 mL	0.2026 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481