


## 9-epi-6-Methoxygeniposidic acid

## Chemical Properties

CAS No. :	1215178-87-2
Formula:	C <sub>17</sub> H <sub>24</sub> O <sub>11</sub>
Molecular Weight:	404.36
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	9-epi-6-Methoxygeniposidic acid is a natural product that can be used in related research in the field of life sciences. Its product number is TN8508.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.473 mL	12.3652 mL	24.7304 mL
5 mM	0.4946 mL	2.473 mL	4.9461 mL
10 mM	0.2473 mL	1.2365 mL	2.473 mL
50 mM	0.0495 mL	0.2473 mL	0.4946 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481