


## 5-Acetoxy Ganoderenic acid F

## Chemical Properties

CAS No. : 1309931-91-6  
 Formula: C<sub>32</sub>H<sub>40</sub>O<sub>9</sub>  
 Molecular Weight: 568.65  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	5-Acetoxy Ganoderenic acid F is a natural product that can be used in related research in the field of life sciences. Its product number is TN8527.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7586 mL	8.7928 mL	17.5855 mL
5 mM	0.3517 mL	1.7586 mL	3.5171 mL
10 mM	0.1759 mL	0.8793 mL	1.7586 mL
50 mM	0.0352 mL	0.1759 mL	0.3517 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481