


3'-Methoxycoumestrol

Chemical Properties

CAS No. : 13360-66-2
 Formula: C₁₆H₁₀O₆
 Molecular Weight: 298.25
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	3'-Methoxycoumestrol is a natural product that can be used in related research in the field of life sciences. Its product number is TN8534.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3529 mL	16.7645 mL	33.5289 mL
5 mM	0.6706 mL	3.3529 mL	6.7058 mL
10 mM	0.3353 mL	1.6764 mL	3.3529 mL
50 mM	0.0671 mL	0.3353 mL	0.6706 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481