


## Congmuyenoside I

## Chemical Properties

CAS No. : 1417579-26-0  
 Formula: C<sub>48</sub>H<sub>78</sub>O<sub>19</sub>  
 Molecular Weight: 958.00  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Congmuyenoside I is a natural product that can be used in related research in the field of life sciences. Its product number is TN8555.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0438 mL	5.2192 mL	10.4384 mL
5 mM	0.2088 mL	1.0438 mL	2.0877 mL
10 mM	0.1044 mL	0.5219 mL	1.0438 mL
50 mM	0.0209 mL	0.1044 mL	0.2088 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481