


Notoginsenoside ST4

Chemical Properties

CAS No. : 155683-02-6
 Formula: C47H80O17
 Molecular Weight: 917.14
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|-------------|--|
| Description | Notoginsenoside ST4 is a natural product that can be used in related research in the field of life sciences. Its product number is TN8574. |
|-------------|--|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.0903 mL | 5.4517 mL | 10.9035 mL |
| 5 mM | 0.2181 mL | 1.0903 mL | 2.1807 mL |
| 10 mM | 0.109 mL | 0.5452 mL | 1.0903 mL |
| 50 mM | 0.0218 mL | 0.109 mL | 0.2181 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481