

Salidroside-3-O-glucopyranoside

Chemical Properties

CAS No. : 1642575-96-9

Formula: C₂₀H₃₀O₁₃


Molecular Weight: 478.45

Keep away from direct sunlight, Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Salidroside-3-O-glucopyranoside is an anti-inflammatory natural product exhibiting inhibitory activity against TNF (IC ₅₀ =38.18 mg/L).
Targets(IC ₅₀)	TNF

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0901 mL	10.4504 mL	20.9008 mL
5 mM	0.418 mL	2.0901 mL	4.1802 mL
10 mM	0.209 mL	1.045 mL	2.0901 mL
50 mM	0.0418 mL	0.209 mL	0.418 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Xu L, Luo M, Lin LM, et al. Three new phenolic glycosides from the Tibetan medicinal plant Aconitum tanguticum. J Asian Nat Prod Res. 2013;15(7):743-749. doi:10.1080/10286020.2013.799145

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481