


Notoginsenoside A

Chemical Properties

CAS No. :	193895-21-5
Formula:	C54H92O24
Molecular Weight:	1125.29
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Notoginsenoside A is a natural product that can be used in related research in the field of life sciences. Its product number is TN8618.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8887 mL	4.4433 mL	8.8866 mL
5 mM	0.1777 mL	0.8887 mL	1.7773 mL
10 mM	0.0889 mL	0.4443 mL	0.8887 mL
50 mM	0.0178 mL	0.0889 mL	0.1777 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481