


## Delphinidin 3-O-galactoside

## Chemical Properties

CAS No. :	197250-28-5
Formula:	C <sub>21</sub> H <sub>21</sub> O <sub>12</sub>
Molecular Weight:	465.39
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Delphinidin 3-O-galactoside is a natural product that can be used in related research in the field of life sciences. Its product number is TN8622.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1487 mL	10.7437 mL	21.4874 mL
5 mM	0.4297 mL	2.1487 mL	4.2975 mL
10 mM	0.2149 mL	1.0744 mL	2.1487 mL
50 mM	0.043 mL	0.2149 mL	0.4297 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481