


Ethyl ganoderenate D

Chemical Properties

CAS No. :	2226858-23-5
Formula:	C ₃₂ H ₄₄ O ₇
Molecular Weight:	540.70
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ethyl ganoderenate D is a natural product that can be used in related research in the field of life sciences. Its product number is TN8641.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8495 mL	9.2473 mL	18.4945 mL
5 mM	0.3699 mL	1.8495 mL	3.6989 mL
10 mM	0.1849 mL	0.9247 mL	1.8495 mL
50 mM	0.037 mL	0.1849 mL	0.3699 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481