

## 1,6,8-trihydroxy-2,7-dimethoxy-3- methylanthraquinone

## Chemical Properties


CAS No. : 2366153-27-5

Formula: C17H14O7

Molecular Weight: 330.29

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

1,6,8-trihydroxy-2,7-dimethoxy-3- methylanthraquinone is a natural product that can be used in related research in the field of life sciences. Its product number is TN8649.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0276 mL	15.1382 mL	30.2764 mL
5 mM	0.6055 mL	3.0276 mL	6.0553 mL
10 mM	0.3028 mL	1.5138 mL	3.0276 mL
50 mM	0.0606 mL	0.3028 mL	0.6055 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481