


Qinyangshengenin-3-O-glycoside

Chemical Properties

CAS No. :	250217-53-9
Formula:	C34H46O11
Molecular Weight:	630.72
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Qinyangshengenin-3-O-glycoside is a natural product that can be used in related research in the field of life sciences. Its product number is TN8657.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5855 mL	7.9274 mL	15.8549 mL
5 mM	0.3171 mL	1.5855 mL	3.171 mL
10 mM	0.1585 mL	0.7927 mL	1.5855 mL
50 mM	0.0317 mL	0.1585 mL	0.3171 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481