


## Petunidin 3,5-glucoside

## Chemical Properties

CAS No. : 25846-73-5  
 Formula: C<sub>28</sub>H<sub>33</sub>ClO<sub>17</sub>  
 Molecular Weight: 677.00  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Petunidin 3,5-glucoside is a natural product that can be used in related research in the field of life sciences. Its product number is TN8660.
-------------	------------------------------------------------------------------------------------------------------------------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4771 mL	7.3855 mL	14.771 mL
5 mM	0.2954 mL	1.4771 mL	2.9542 mL
10 mM	0.1477 mL	0.7386 mL	1.4771 mL
50 mM	0.0295 mL	0.1477 mL	0.2954 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481