


## Gyenoside XCI

## Chemical Properties

CAS No. : 2793436-99-2  
 Formula: C47H80O17  
 Molecular Weight: 916.00  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Gyenoside XCI is a natural product that can be used in related research in the field of life sciences. Its product number is TN8665.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0917 mL	5.4585 mL	10.917 mL
5 mM	0.2183 mL	1.0917 mL	2.1834 mL
10 mM	0.1092 mL	0.5459 mL	1.0917 mL
50 mM	0.0218 mL	0.1092 mL	0.2183 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481