

3'-O-Methy Isosaponarin

Chemical Properties


CAS No. : 2803327-85-5

Formula: C₂₈H₃₂O₁₆

Molecular Weight: 624.50

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

3'-O-Methy Isosaponarin is a natural product that can be used in related research in the field of life sciences. Its product number is TN8666.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6013 mL	8.0064 mL	16.0128 mL
5 mM	0.3203 mL	1.6013 mL	3.2026 mL
10 mM	0.1601 mL	0.8006 mL	1.6013 mL
50 mM	0.032 mL	0.1601 mL	0.3203 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481