


24(28)-Dehydroergosterol

Chemical Properties

CAS No. : 29560-24-5
 Formula: C₂₈H₄₂O
 Molecular Weight: 394.64
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description 24(28)-Dehydroergosterol is a natural product that can be used in related research in the field of life sciences. Its product number is TN8672.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.534 mL	12.6698 mL	25.3395 mL
5 mM	0.5068 mL	2.534 mL	5.0679 mL
10 mM	0.2534 mL	1.267 mL	2.534 mL
50 mM	0.0507 mL	0.2534 mL	0.5068 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481