


## Telephiose A

## Chemical Properties

CAS No. :	297749-29-2
Formula:	C <sub>36</sub> H <sub>44</sub> O <sub>20</sub>
Molecular Weight:	796.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Telephiose A is a natural product that can be used in related research in the field of life sciences. Its product number is TN8673.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2559 mL	6.2795 mL	12.559 mL
5 mM	0.2512 mL	1.2559 mL	2.5118 mL
10 mM	0.1256 mL	0.628 mL	1.2559 mL
50 mM	0.0251 mL	0.1256 mL	0.2512 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481