


## Telephiose C

## Chemical Properties

CAS No. : 297749-31-6  
 Formula: C32H40O18  
 Molecular Weight: 712.22  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	Telephiose C is a natural product that can be used in related research in the field of life sciences. Its product number is TN8674.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4041 mL	7.0203 mL	14.0406 mL
5 mM	0.2808 mL	1.4041 mL	2.8081 mL
10 mM	0.1404 mL	0.702 mL	1.4041 mL
50 mM	0.0281 mL	0.1404 mL	0.2808 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481