


## Hesperetin Dihydrochalcone

## Chemical Properties

CAS No. :	35400-60-3
Formula:	C <sub>16</sub> H <sub>16</sub> O <sub>6</sub>
Molecular Weight:	304.30
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Hesperetin Dihydrochalcone is a natural product that can be used in related research in the field of life sciences. Its product number is TN8682.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2862 mL	16.4312 mL	32.8623 mL
5 mM	0.6572 mL	3.2862 mL	6.5725 mL
10 mM	0.3286 mL	1.6431 mL	3.2862 mL
50 mM	0.0657 mL	0.3286 mL	0.6572 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481