

Eleutheroside I

Chemical Properties


CAS No. : 35790-94-4

Formula: C41H66O11

Molecular Weight: 734.97

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Eleutheroside I is a natural product that can be used in related research in the field of life sciences. Its product number is TN8684.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3606 mL	6.803 mL	13.606 mL
5 mM	0.2721 mL	1.3606 mL	2.7212 mL
10 mM	0.1361 mL	0.6803 mL	1.3606 mL
50 mM	0.0272 mL	0.1361 mL	0.2721 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481