


Ergosta-7,22,24(28)-trien-3 β -ol

Chemical Properties

CAS No. : 36680-39-4
 Formula: C₂₈H₄₄O
 Molecular Weight: 396.66
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ergosta-7,22,24(28)-trien-3 β -ol is a natural product that can be used in related research in the field of life sciences. Its product number is TN8686.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5211 mL	12.6053 mL	25.2105 mL
5 mM	0.5042 mL	2.5211 mL	5.0421 mL
10 mM	0.2521 mL	1.2605 mL	2.5211 mL
50 mM	0.0504 mL	0.2521 mL	0.5042 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481