


## 1-erucoyl-sn-glycero-3-phosphocholine

## Chemical Properties

CAS No. : 503271-85-0  
 Formula: C30H60NO7P  
 Molecular Weight: 577.78  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-erucoyl-sn-glycero-3-phosphocholine is a natural product that can be used in related research in the field of life sciences. Its product number is TN8703.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7308 mL	8.6538 mL	17.3076 mL
5 mM	0.3462 mL	1.7308 mL	3.4615 mL
10 mM	0.1731 mL	0.8654 mL	1.7308 mL
50 mM	0.0346 mL	0.1731 mL	0.3462 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481