


## 4,3',5'-Trihydroxyresveratrol

## Chemical Properties

CAS No. :	637776-83-1
Formula:	C <sub>14</sub> H <sub>12</sub> O <sub>6</sub>
Molecular Weight:	276.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4,3',5'-Trihydroxyresveratrol is a natural product that can be used in related research in the field of life sciences. Its product number is TN8733.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.620 mL	18.1002 mL	36.2004 mL
5 mM	0.724 mL	3.620 mL	7.2401 mL
10 mM	0.362 mL	1.810 mL	3.620 mL
50 mM	0.0724 mL	0.362 mL	0.724 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481