


## Hemerocallal A

## Chemical Properties

|                   |   |
|-------------------|---|
| CAS No. :         | 659747-79-2   |
| Formula:          | C <sub>20</sub> H <sub>28</sub> O <sub>2</sub>  |
| Molecular Weight: | 300.44  |
| Storage:          | Powder: -20°C for 3 years   In solvent: -80°C for 1 year<br>Actual storage temperature shall be subject to the COA. |



## Biological Description

|             |   |
|-------------|---|
| Description | Hemerocallal A is a natural product that can be used in related research in the field of life sciences. Its product number is TN8742. |
|-------------|---|

## Preparing Stock Solutions

|       | 1mg       | 5mg        | 10mg       |
|-------|-----------|------------|------------|
| 1 mM  | 3.3285 mL | 16.6423 mL | 33.2845 mL |
| 5 mM  | 0.6657 mL | 3.3285 mL  | 6.6569 mL  |
| 10 mM | 0.3328 mL | 1.6642 mL  | 3.3285 mL  |
| 50 mM | 0.0666 mL | 0.3328 mL  | 0.6657 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481