

Veratrine

Chemical Properties

CAS No. :	8051-02-3
Formula:	C ₃₂ H ₄₉ N _O ₉
Molecular Weight:	591.74
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Veratrine is a natural product that can be used in related research in the field of life sciences. Its product number is TN8758.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6899 mL	8.4497 mL	16.8993 mL
5 mM	0.338 mL	1.6899 mL	3.3799 mL
10 mM	0.169 mL	0.845 mL	1.6899 mL
50 mM	0.0338 mL	0.169 mL	0.338 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481