


Salannal

Chemical Properties

CAS No. :	86160-86-3
Formula:	C34H44O10
Molecular Weight:	612.71
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Salannal is a natural product that can be used in related research in the field of life sciences. Its product number is TN8768.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6321 mL	8.1605 mL	16.3209 mL
5 mM	0.3264 mL	1.6321 mL	3.2642 mL
10 mM	0.1632 mL	0.816 mL	1.6321 mL
50 mM	0.0326 mL	0.1632 mL	0.3264 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481