


Gyenoside H-2

Chemical Properties

CAS No. : 862286-47-3
 Formula: C₄₉H₇₈O₁₈
 Molecular Weight: 955.15
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Gyenoside H-2 is a natural product that can be used in related research in the field of life sciences. Its product number is TN8770.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.047 mL	5.2348 mL	10.4696 mL
5 mM	0.2094 mL	1.047 mL	2.0939 mL
10 mM	0.1047 mL	0.5235 mL	1.047 mL
50 mM	0.0209 mL	0.1047 mL	0.2094 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481