


Herpetotriol

Chemical Properties

CAS No. :	89020-08-6
Formula:	C30H32O9
Molecular Weight:	536.57
Storage:	Keep away from direct sunlight, Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Herpetotriol is a natural lignan with potential hepatoprotective effects which can be employed in studies of jaundiced hepatitis and dyspepsia.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8637 mL	9.3184 mL	18.6369 mL
5 mM	0.3727 mL	1.8637 mL	3.7274 mL
10 mM	0.1864 mL	0.9318 mL	1.8637 mL
50 mM	0.0373 mL	0.1864 mL	0.3727 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481