


Tetraacetylisorhamnetin

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 990-19-2 |
| Formula: | C ₂₄ H ₂₀ O ₁₁ |
| Molecular Weight: | 484.41 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|-------------|--|
| Description | Tetraacetylisorhamnetin is a natural product that can be used in related research in the field of life sciences. Its product number is TN8805. |
|-------------|--|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.0644 mL | 10.3218 mL | 20.6437 mL |
| 5 mM | 0.4129 mL | 2.0644 mL | 4.1287 mL |
| 10 mM | 0.2064 mL | 1.0322 mL | 2.0644 mL |
| 50 mM | 0.0413 mL | 0.2064 mL | 0.4129 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481