

Moldapyranone B

Chemical Properties

CAS No. :

Formula: C₁₄H₁₂O₅

Molecular Weight: 260.25

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Moldapyranone B is a natural product that can be used in related research in the field of life sciences. Its product number is TN8809.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8425 mL	19.2123 mL	38.4246 mL
5 mM	0.7685 mL	3.8425 mL	7.6849 mL
10 mM	0.3842 mL	1.9212 mL	3.8425 mL
50 mM	0.0768 mL	0.3842 mL	0.7685 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481