

Naringin chalcone trimethyl ether

Chemical Properties

CAS No. :

Formula: C₃₀H₃₈O₁₄

Molecular Weight: 622.62

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Naringin chalcone trimethyl ether is a natural product that can be used in related research in the field of life sciences. Its product number is TN8815.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6061 mL	8.0306 mL	16.0612 mL
5 mM	0.3212 mL	1.6061 mL	3.2122 mL
10 mM	0.1606 mL	0.8031 mL	1.6061 mL
50 mM	0.0321 mL	0.1606 mL	0.3212 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481