

Altissimacoumarin K

Chemical Properties

CAS No. :	2555262-76-3
Formula:	C ₂₁ H ₂₈ O ₈
Molecular Weight:	408.44
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Altissimacoumarin K is a natural product that can be used in related research in the field of life sciences. Its product number is TN8826.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4483 mL	12.2417 mL	24.4834 mL
5 mM	0.4897 mL	2.4483 mL	4.8967 mL
10 mM	0.2448 mL	1.2242 mL	2.4483 mL
50 mM	0.049 mL	0.2448 mL	0.4897 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481