

## Dihydropyrazine

## Chemical Properties

CAS No. :

Formula: C<sub>4</sub>H<sub>6</sub>N<sub>2</sub>

Molecular Weight: 82.10

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Dihydropyrazine is a natural product that can be used in related research in the field of life sciences. Its product number is TN8839.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	12.1803 mL	60.9013 mL	121.8027 mL
5 mM	2.4361 mL	12.1803 mL	24.3605 mL
10 mM	1.218 mL	6.0901 mL	12.1803 mL
50 mM	0.2436 mL	1.218 mL	2.4361 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481