

Bacopaside VI

Chemical Properties

CAS No. :

Formula: C41H66O16S

Molecular Weight: 847.00

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bacopaside VI is a natural product that can be used in related research in the field of life sciences. Its product number is TN8852.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1806 mL	5.9032 mL	11.8064 mL
5 mM	0.2361 mL	1.1806 mL	2.3613 mL
10 mM	0.1181 mL	0.5903 mL	1.1806 mL
50 mM	0.0236 mL	0.1181 mL	0.2361 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481