

Human IgG1(E356D/M358L/L234A/L235A)-kappa Isotype control

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Human IgG1(E356D/M358L/L234A/L235A)-kappa Isotype control is the isotype control for Human IgG1(E356D/M358L/L234A/L235A)-kappa antibody.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481