

Human IgG4(S228P/D265A/L309V/R409K)-kappa Isotype control

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Human IgG4(S228P/D265A/L309V/R409K)-kappa Isotype control is the isotype control for Human IgG4(S228P/D265A/L309V/R409K)-kappa antibody.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481