

## Human IgG1(E356D/M358L/L234A/L235A/F405L)-kappa Isotype control

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Human IgG1(E356D/M358L/L234A/L235A/F405L)-kappa Isotype control is the isotype control for Human IgG1(E356D/M358L/L234A/L235A/F405L)-kappa antibody.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481