

Mouse Anti-Human IgA Secondary Antibody-HRP (7L75)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	HRP
Clone:	7L75
Purification:	Protein A

Applications

Verified Activity:	ELISA of human immunoglobulins shows Mouse Anti-Human IgA Secondary Antibody (HRP) reacts to human IgA, recombinant human IgA1. No cross reactivity with human IgG, IgE. The plate was coated with 1 µg/mL of different immunoglobulins. 0.5 µg/mL of Mouse Anti-Human IgA Secondary Antibody (HRP) was used as the secondary antibody.
Application:	ELISA
Recommended	ELISA: 0.1-1 µg/mL

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Native human IgA
Antigen Species:	Human
Synonyms:	IgA1;AL928768.3

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481