

## Anti-Enterokinase Antibody-HRP (9C751)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Bovine
Conjugation:	HRP
Clone:	9C751
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 0.1-1 µg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Bovine Enterokinase / PRSS7 Protein (TMPY-00274)
Antigen Species:	Bovine
Synonyms:	transmembrane protease, serine 15
Biology Area:	Serine Proteases and Regulators

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481