

Anti-Glypican 3/GPC3 Antibody-PE (4J109)

Product Details

| | |
|---------------|------------|
| Ig Type: | Rabbit IgG |
| Reactivity: | Human |
| Conjugation: | PE |
| Clone: | 4J109 |
| Purification: | Protein A |

Applications

| | |
|--------------------|--|
| Verified Activity: | Flow cytometric analysis of Human GPC3 expression on HepG2 cells. Cells were stained with PE-conjugated anti-Human GPC3. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells. |
| Application: | FCM |
| Recommended | 10 µl/Test, 0.1 mg/ml |

Properties

| | |
|----------------------|--|
| Stability & Storage: | Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|---|
| Immunogen: | A synthetic peptide: center region of the Human Glypican 3 / GPC3 / OCI-5 |
| Antigen Species: | Human |
| Synonyms: | SGBS1;MXR7;SDYS;DGSX;SGB;Glypican 3;SGBS;GTR2-2;OCI-5 |
| Biology Area: | Cancer Drug Targets |

Research Background

Glypican-3, also known as Intestinal protein OCI-5, GPC3, and OCI5, is a member of the glypican family. It belongs to the glypican family and is highly expressed in the lung, liver, and kidney. It is a heparan sulfate proteoglycan, which is overexpressed in various neoplasms such as hepatocellular carcinoma, malignant melanoma, and testicular yolk sac tumor, and plays an important role in cell growth and differentiation. GPC3 function is tissue-dependent. In some tissues, GPC3 acts as a tumor suppressor gene, whereas in others, it acts as an oncofetal protein. Studies have shown that GPC3 is a reliable marker for hepatocellular carcinoma. The sensitivity and specificity exceed both alpha-fetoprotein and hepatocyte-paraffin1. GPC3 immunohistochemistry can aid in the differentiation of testicular germ cell tumors, being expressed in all yolk sac tumors but not in seminomas. GPC3 expression has also been identified in some squamous cell carcinomas of the lung and clear cell carcinomas of the ovary. The role of GPC3 in melanomas is still controversial. Thus, Glypican-3 is currently regarded as a tumor marker and potential target for immunotherapy. Cancer Immunotherapy Immune Checkpoint Immunotherapy Targeted Therapy

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481