

Anti-CCR1/CD191 Antibody (1E968)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1E968
Purification:	Protein A

Applications

Verified Activity:	<ol style="list-style-type: none">1. Immunochemical staining of human CCR1 in human spleen with rabbit monoclonal antibody (1:500, formalin-fixed paraffin embedded sections).2. Immunochemical staining of human CCR1 in human tonsil with rabbit monoclonal antibody (1:500, formalin-fixed paraffin embedded sections).
Application:	IHC-P
Recommended	IHC-P: 1:250-1:1000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: center region of the Human CCR1
Antigen Species:	Human
Synonyms:	CCR-1;CKR-1;MIP1aR;CKR1;CMKBR1;CMKR1;CD191;HM145;SCYAR1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481