

## Anti-CCR2/CD192 Antibody (30968)

## Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	30968
Purification:	Protein A

## Applications

Verified Activity:	Flow cytometric analysis of Human CD192 expression on human whole blood monocytes. Cells were stained with purified anti-Human CD192, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.
Application:	FCM
Recommended	FCM: 1:25-1:100

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	A synthetic peptide: center region of the Human CCR2/CD192
Antigen Species:	Human
Synonyms:	CCR2A;CC-CKR-2;MCP-1-R;CKR2A;CMKBR2;CKR2;CKR2B;CD192;CCR-2;PCLUD;CCR2B;FLJ78302
Biology Area:	Neuroinflammation

## Research Background

The C-C motif chemokine receptor-2 (CCR2) is a G protein-coupled receptor (GPCR), made up of a carboxy-terminus, extracellular amino terminus, and a hydrophobic transmembrane domain consisting of 7 amino acid segments that mediates agonist-dependent calcium mobilization and inhibition of adenylyl cyclase. CCR2 is expressed on monocytes, immature dendritic cells, and T-cell subpopulations, and mediates their migration towards endogenous CCL2. CCR2 is necessary for macrophage-dependent inflammatory responses and the development of atherosclerosis. In mice, CCR2 deficiency reduced macrophage content, increased adiponectin expression, ameliorated hepatic steatosis, and improved systemic glucose homeostasis and insulin sensitivity. Resistance to HIV-1 infection or delayed progression to AIDS may be linked to CCR2 polymorphisms. Furthermore, CCR2 mRNA was highly expressed in prostate cancer (PCa) metastatic tissues compared with benign prostate tissues according to real-time RT-PCR3, suggesting that CCR2 may contribute to PCa development.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481