

## Anti-IL-27 Antibody (8P954)

## Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8P954
Purification:	Protein A

## Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Human IL-27 Protein (TMPY-03142)
Antigen Species:	Human
Synonyms:	IL-27;interleukin 27;p28;IL30

## Research Background

IL-27 protein is a member of the IL-6 superfamily, which is expressed on monocytes, endothelial cells, and dendritic cells. IL-27 protein is also referred to as the IL-12 p35-related protein, p28, is one subunit of a heterodimeric cytokine complex, and associates with another subunit EB13 (EBV-induced gene 3), and IL-12 p40-related protein (IL-27B). IL-27 protein is an early product of activated antigen-presenting cells and drives the rapid clonal expansion of naive CD4(+) T cells and plays a role in the early regulation of Th1 cells initiation which drives efficient adaptive immune response. IL-27 protein has an antiproliferative activity on melanomas through WSX-1/STAT1 signaling. Thus, IL-27 protein may be an attractive candidate as an antitumor agent applicable to cancer immunotherapy.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481