

Anti-IL-1F10 Antibody (6N340)

Product Details

Ig Type:	Mouse IgG2b
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6N340
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human IL1F10 / IL-38 Protein (TMPY-04175)
Antigen Species:	Human
Synonyms:	FKSG75;IL1-theta;IL1-θ;FIL1-theta;IL-38;interleukin 1 family, member 10 (theta);IL-1HY2;IL1HY2;interleukin 1 family, member 10 (θ);FIL1-θ

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481