

Anti-IL-36 gamma Antibody (7U281)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7U281
Purification:	Protein A

Applications

Verified Activity:	<ol style="list-style-type: none">1. Immunochemical staining of human IL1F9 in human lung with rabbit monoclonal antibody (1:20000, formalin-fixed paraffin embedded sections).2. Immunochemical staining of human IL1F9 in human esophagus with rabbit monoclonal antibody (1:20000, formalin-fixed paraffin embedded sections).
Application:	IHC-P
Recommended	IHC-P: 1:10000-1:20000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human IL36G / IL1F9 Protein (TMPY-02852)
Antigen Species:	Human
Synonyms:	interleukin 36, γ ;interleukin 36, gamma;IL36 γ ;IL36g;Il1f9

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481