

## Anti-IGF1R/CD221 Antibody-FITC (1V195)

## Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	FITC
Clone:	1V195
Purification:	Protein A

## Applications

Verified Activity:	Flow cytometric analysis of anti-Human IGF1R on HeLa cells.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

## Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Human IGF1R / CD221 protein (TMPY-00834)
Antigen Species:	Human
Synonyms:	IGF-I R;insulin-like growth factor 1 receptor;IGFR;JTK13;IGFIR;CD221;IGF1 Receptor
Biology Area:	Receptor Tyrosine Kinases (RTKs), Cancer Drug Targets

## Research Background

The insulin-like growth factor-1 receptor (IGF1R) is a transmembrane tyrosine kinase involved in several biological processes including cell proliferation, differentiation, DNA repair, and cell survival. This a disulfide-linked heterotetrameric transmembrane protein consisting of two  $\alpha$  and two  $\beta$  subunits, and among which, the  $\alpha$  subunit is extracellular while the  $\beta$  subunit has an extracellular domain, a transmembrane domain, and a cytoplasmic tyrosine kinase domain. The IGF1R signaling pathway is activated in the mammalian nervous system from the early developmental stages. Its major effect on developing neural cells is to promote their growth and survival. This pathway can integrate its action with signaling pathways of growth and morphogenetic factors that induce cell fate specification and selective expansion of specified neural cell subsets. Modulation of cell migration is another possible role that IGF1R activation may play in neurogenesis. In the mature brain, IGF-I binding sites have been found in different regions of the brain, and multiple reports confirmed a strong neuroprotective action of the IGF-IR against different pro-apoptotic insults. IGF1R is an important signaling molecule in cancer cells and plays an essential role in the establishment and maintenance of the transformed phenotype. Inhibition of IGF1R signaling thus appears to be a promising strategy to interfere with the growth and survival of cancer cells. IGF1R is frequently overexpressed by tumors and mediates proliferation and apoptosis protection. IGF signaling also influences hypoxia signaling, protease secretion, tumor cell motility, and adhesion, and thus can affect the propensity for invasion and

## A DRUG SCREENING EXPERT

---

metastasis. Therefore, IGF1R is now an attractive anti-cancer treatment target. Cancer Immunotherapy Immune Checkpoint Immunotherapy Targeted Therapy

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481