

## Anti-Galectin-1 Antibody (5G581)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5G581
Purification:	Protein A

### Applications

Verified Activity:	<p>1. Anti-Galectin1 mouse monoclonal antibody at 1:500</p> <ul style="list-style-type: none"><li>-Lane A: Hela Whole Cell Lysate.</li><li>-Lane B: A431 Whole Cell lysate.</li><li>-Lysates/proteins at 30 µg per lane.</li><li>-Secondary</li><li>-Goat Anti-Mouse IgG H&amp;L (Dylight800) at 1/15000 dilution.</li><li>-Developed using the Odyssey technique.</li><li>-Performed under reducing conditions.</li></ul> <p>2. Galectin1 was immunoprecipitated using:</p> <ul style="list-style-type: none"><li>-Lane A:0.5 mg Hela Whole Cell Lysate</li><li>-0.5 µL anti-Galectin1 mouse monoclonal antibody and 15 µl of 50 % Protein G agarose.</li><li>-Primary antibody:</li><li>-Anti-Galectin1 mouse monoclonal antibody, at 1:250 dilution.</li><li>-Secondary antibody:</li><li>-Dylight 800-labeled antibody to Mouse IgG (H+L), at 1:7500 dilution.</li><li>-Developed using the odyssey technique.</li><li>-Performed under reducing conditions.</li><li>-Predicted band size: 15 kDa.</li><li>-Observed band size: 15 kDa</li></ul>
Application:	IP,WB
Recommended	WB: 1:500-1:1000; IP: 0.2-1 µL/mg of lysate

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen: Recombinant Protein: Human Galectin-1 / LGALS1 protein (TMPY-00976)

Antigen Species: Human

Synonyms: GAL1;GBP;lectin, galactoside-binding, soluble, 1;Galectin-1

---

### Research Background

Galectin-1 (Gal-1, GAL1), is a member of the galectins, a family of animal lectins ranging from *Caenorhabditis elegans* to humans, which is defined by their affinity for beta-galactosides and by significant sequence similarity in the carbohydrate-binding site. It is a homodimer with a subunit molecular mass of 14.5 kDa, which contains six cysteine residues per subunit. The cysteine residues should be in a free state to maintain a molecular structure that is capable of showing lectin activity. This endogenous lectin widely expressed at sites of inflammation and tumor growth has been postulated as an attractive immunosuppressive agent to restore immune cell tolerance and homeostasis in autoimmune and inflammatory settings. On the other hand, galectin-1 contributes to different steps of tumor progression including cell adhesion, migration, and tumor-immune escape, suggesting that blockade of galectin-1 might result in therapeutic benefits in cancer. Several potential glycoprotein ligands for galectin-1 have been identified, including lysosome-associated membrane glycoproteins and fibronectin, laminin, as well as T-cell glycoproteins CD43 and CD45. Evidence points to Gal-1 and its ligands as one of the master regulators of such immune responses as T-cell homeostasis and survival, T-cell immune disorders, inflammation, and allergies as well as host-pathogen interactions.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481