

Anti-Serpin A3 Antibody (9Y470)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9Y470
Purification:	Protein A

Applications

Application:	ELISA,ELISA(Det)
Recommended	ELISA: 1:5000-1:10000; ELISA(Det): 1:1000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human SerpinA3 Protein (TMPY-01303)
Antigen Species:	Human
Synonyms:	serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 3;SERPINA3; serpin peptidase inhibitor, clade A (α -1 antiproteinase, antitrypsin), member 3;GIG25; MGC88254;AACT;ACT;GIG24

Research Background

SerpinA3, also known as Alpha 1-antichymotrypsin (AACT), is a plasma alpha globulin glycoprotein, and is a member of serpin superfamily of the serine protease inhibitors consisting of at least 35 members. SerpinA3 has been demonstrated to inhibit the activity of certain serine proteases, such as cathepsin G found in neutrophils, and chymases present in mast cells, by inducing a major conformational rearrangement, and thus protects some tissues from damage caused by proteolytic enzymes. This enzyme is produced primarily in the liver, and is identified as an acute-phase inflammatory protein. SerpinA3 deficiency has been associated with liver disease, and mutations of this gene have been observed in patients with Parkinson disease and chronic obstructive pulmonary disease. Besides, ACT gene polymorphism has been implicated with Alzheimer's disease (AD), cerebral amyloid angiopathy (CAA), as well as stroke, since SerpinA3 is a major constituent of the plaques in AD and an inhibitor of amyloid beta peptide degradation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481